

Disposable Digital Forehead Thermometer Products

Instructions for Use

How to Read

1. Remove backing liner. Apply to center of forehead, just above the eyebrows. Read the indicator when color stops changing and while it is still on the forehead.
2. Green indicates the correct temperature; if green does not appear, the temperature is between blue and tan.

Intended Use

- Measure human forehead surface temperature and monitor temperature trends

Precautions

- Environmental temperatures might affect the performance of this product
- Discard after use

Liquid Crystal thermometers are to be used on the forehead area only and estimate the core body temperature by measuring the temperature of the forehead. It is an ideal way to check for changes in temperature. A Liquid Crystal Forehead Thermometer is not as accurate for determining the actual core body temperature as a clinical (oral / rectal) thermometer. Instead it is designed to be an unbreakable, quick reading device for temperature screening between uses of a clinical thermometer.

Forehead skin temperatures are lower than oral temperatures.

Accurate to $\pm 1.0^{\circ}\text{F}$ / $\pm 0.6^{\circ}\text{C}$, the readings on a Liquid Crystal Forehead Thermometer have been adjusted to approximate oral temperature readings during normal to feverish conditions. If you think your temperature is different from that shown by a Liquid Crystal Forehead Thermometer, you should check by using a clinical thermometer. If fever is present, consult your physician for appropriate treatment.

For External use only

This product contains no latex

Storage / Handling:

Products should be stored out of direct sunlight, dampness, or heat

Cleaning:

Wipe clean using soap and water or rubbing alcohol. Do not soak.

Conforms to ASTM E-1061-01



LCR Hallcrest, 1911 Pickwick Lane, Glenview IL 60026 USA



LCR Hallcrest Ltd. Dock Road, Connah's Quay, UK CH5 4DS

Notified Body: Det Norske Veritas 0434